

# FIBER-OPTICAL SERIES

USB

- > WIDE DYNAMIC RANGE
- > WIDE EXPOSURE RANGE



**HDR** - 500/600

Digital Dental X-ray Imaging System

CE 0123

## F Functional Features

### Fiber Optic Plate

The built-in FOP can take more realistic high-definition images, and even cracked teeth can be easily found.



HDR-500

lower molars




HDR-600

lower molars

### Wide Exposure Range

Size 1.3 sensor, shooting width 22.5mm



HDR-500

lower front teeth



HDR-600


lower front teeth

### Wide Dynamic Range

Easy to shoot from low dose to high dose, greatly reducing the requirements for filming experience. Image clarity and sensitivity are significantly improved.

\*Pictures from the imaging output

## P Product Parameters



Shanghai Handy Medical Equipment Co., Ltd  
Tel: +86 021-36161602 / +86 021-56766108  
✉ Business Contact: sale@handycreate.com

Website: www.handyimaging.com  
✉ After Sales Service: support@handycreate.com