



GOOD SMILE DESIGN

PORTABLE

DIGITAL IMAGING PLATE SCANNER

COMPACT SIZE WITH POWERFUL FEATURES



4 Sizes (0/1/2/3) Imaging Plates Available



1.5kg Lightweight



Imaging Within 6s

HDS -500

Digital Imaging Plate Scanner



F Functional Features

Mini and Portable

With weight < 1.5kg, it is highly integrated, ultra small and convenient to be moved.

One-Click Imaging


It can output images within 6s, meeting your fast imaging requirement.

Flexible Imaging Plates

It is easier to fit the teeth, and the shooting is more comfortable and accurate.

4 Sizes

It is suitable for size 0/1/2/3 imaging plates to meet different application scenarios.



HDS-500 Size 2 *Pictures from clinic application




HDS-500 Size 2 *Pictures from clinic application

P Product Parameters

Laser Spot Size	Imaging Time	Laser Wavelength	Weight	ADC	Operating System
35μm	≤6sec	660nm	<1.5kgs	14bit	win7/10 (32bit&64bit)

S Size Show



Shanghai Handy Medical Equipment Co., Ltd

Tel: +86 021-36161602 / +86 021-56766108

Business Contact: sale@handycreate.com

Website: www.handyimaging.com

After Sales Service: support@handycreate.com